# SOLE™

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**Expect an adjustment period.** It is normal for some people to experience discomfort as your feet adjust to the sensation of our arch support. Start by gradually introducing the footbeds over a few days, or heat mold if necessary. Stick with us, and soon your feet will feel better than ever.

### Follow these steps for ultimate comfort:

**1 Remove the factory insole from your footwear if possible.**

**2 Trim your SOLE footbeds if needed.**

Trace the factory insole from your footwear as a guide.

Only trim the forefoot area, and don't trim the footbed too short!

**3 Custom mold your footbeds.**

#### Wear Mold

Slowly break-in your footbeds over 3 to 7 days. Gradually increase the duration you wear them for each day, and avoid high-intensity activities. If you find the arch too firm, follow the heat molding instructions to ease the transition and speed up a comfortable fit.

#### Heat Mold

**Always allow the oven to fully preheat before inserting footbeds.** Pay close attention and only heat footbeds for the stipulated time at the stipulated temperature. Overheating your footbeds can cause melting or cracking.

**1 Pre-heat your oven.**

**Conventional oven:**

set to bake at 200°F / 90°C.

**Fan-assisted oven:**

start at 175°F / 80°C.

**Gas oven:**

start at lowest gas mark.

**2 Bake for 2 minutes.**

Place the footbeds on a baking sheet in the pre-heated oven.

Bake for 2 minutes or until the opti-therm indicator turns black.

If the opti-therm indicator has not turned black after 3 minutes, remove the footbeds.

**3 Place footbeds in your shoes.**

Make sure the footbeds are laying flat in your footwear and have not bunched up.

**4 Put your shoes on immediately.**

Make sure your heels are against the back of the shoe.

Lace your footwear up as you normally would.

**5 Stand up straight and stand still for 2 minutes.**

After 2 minutes, your footbeds will be molded to your feet, providing customized support and comfort that's unique to you.

Questions?

We can help.

[customerservice@yoursole.com](mailto:customerservice@yoursole.com)

**Do not use a microwave or air fryer.**

**Health warning:** Diabetics and people with poor circulation or lack of sensation in their feet should consult their health care professional before using this product.

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