# SOLE™

<b>EN</b>	<b>FITTING INSTRUCTIONS</b>	<b>2</b>
<b>FR</b>	<b>INSTRUCTIONS D'AJUSTAGE</b>	<b>3</b>
<b>ES</b>	<b>INSTRUCCIONES DE ADAPTACIÓN</b>	<b>4</b>
<b>DE</b>	<b>ANPASSUNGSANLEITUNG FÜR FUSSBETTEN</b>	<b>5</b>
<b>NO</b>	<b>INNLEGGÅLE TILPASSNINGS INSTRUKSJONER</b>	<b>6</b>
<b>SE</b>	<b>INPASSNING AV FOTBÄDDEN</b>	<b>7</b>

## WEAR MOLD

Remove the factory sock liners from your footwear, insert your new SOLE footbeds into the shoes and wear as normal. The footbed will adapt to your foot after a few days of regular use, giving you a customized fit. If you have low or collapsed arches, we recommend heat molding your footbeds.

## HEAT MOLD

Heat molding will provide the highest level of customization in the shortest amount of time. To heat mold your footbeds, follow these instructions.

- 1 Pre-heat your oven. Do not use a microwave.**  
**Conventional oven:** set to bake at 200°F / 90°C. **Fan-assisted oven:** start at 175°F / 80°C.  
**Gas oven:** start at lowest gas mark.
- 2 If possible, remove the factory-supplied sock liner from your footwear.**  
If the factory-supplied sock liner cannot be removed because it is permanently glued in, leave it in place.
- 3 If necessary, trim your new footbeds to fit.**  
Before you put your new footbeds in the oven, try them in your footwear to make sure they fit properly.
- 4 Bake for two minutes.**  
Place the footbeds in the oven on a baking sheet. Bake for two minutes or until the opti-therm indicator turns black (no more than three minutes).
- 5 Remove from the oven and immediately place inside your footwear.**  
Make sure the footbeds are laying flat in your footwear and have not bunched up.
- 6 Step into your footwear right away.**  
Make sure your heels are against the back of the shoe. Lace your footwear up as you normally would.
- 7 Stand up straight and hold your position for two minutes.**  
After two minutes, your footbeds will be molded to your feet, providing customized support and comfort that's unique to you.

### HEALTH WARNING

Diabetics and people with poor circulation or lack of sensation in their feet should consult their health care professional before using this product.

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